


REVOLUTIONIZE VETERINARY
DENTAL IMAGING WITH

QUICKRAY V-4



Our
Large Format

Veterinary X-ray Sensor



This expanded size allows for a wider field of view, making it easier to capture more anatomy in a single image with sharp, detailed clarity.

Engineered for versatility, this sensor is ideal for a variety of animal sizes and breeds, enabling faster, more efficient diagnostics. It's the upgrade veterinary professionals have been waiting for—powerful, precise, and ready for today's demanding workflows.

Unmatched Image Quality, Unbelievable Value

Your Practice Deserves to Purchase the Right Sensor!

Let our Team of Experts Help You Make the Best Choice

TECHNICAL SPECIFICATIONS

Sensor Tech	CMOS
Dimension	77.1 x 53.8mm
Active Area	67.3 x 46.7 mm
Spatial Resolution	10-20 lp/mm
Pixel Size	18.5 x 18.5um
Pixels	2524x3640 pixels (9.2 Megapixels)
Frame Rate	0.5 fps
Digital Bit Options	14 bit, 16 bit
Connector	USB 2.0
Connectivity	Windows

APPLICATIONS:

- Occlusal
- NDT inspection
- 3D X-ray applications
- General X-ray imaging
- Vet dental

