

QUICKRAY®  
Trusted Since 1998




The  
**HASSLE FREE**  
Intraoral Sensor  
Low Cost, Never Any Hidden Fees

The QuickRay Digital X-Ray Sensor series gives you excellent image clarity at the highest resolution with minimum radiation exposure; and connects directly to your computer's USB port without requiring a special interface box.



TWAIN compatible, QuickRay works in all leading imaging software and can be used in conjunction with existing sensors.

COMPATIBLE WITH



# Your Practice Deserves To Purchase The *Right* Sensor!

Let Our Team of Experts Help You Make the Best Choice

## SENSOR TECHNICAL DATA

|                |  |
|----------------|--|
| SENSOR TECH    | CMOS   |
| ACTIVE AREA    | Size 1: 19.95 x 30.02mm<br>Size 2: 25.99 x 36.00mm |
| DIMENSIONS     | Size 1: 24.4x36.8 mm<br>Size 2: 2 30.5x42.8 mm     |
| BODY THICKNESS | 5.3 mm   |
| PIXEL SIZE     | 19X19 Microns                                      |
| SCINTILLATOR   | CSI  |



## QUICKRAY SENSOR



- Proven Durable Engineering
- Hermetically Sealed Housing
- Shock Absorbing Material
- Ergonomic Design
- Thin with Rounded Corners for Patient Comfort

## CORD



- 2M Cable Length
- Kevlar Reinforced Virtually Unbreakable
- Replaceable
- Direct High Speed 2.0 USB Connectivity



ACTUAL IMAGE

## SENSOR KIT INCLUDES



Protective Sheaths    UNI-VERSE-ALL Positioning Kit    Wall Mount Holder



## Warranty

3 Year Warranty



## Expedited Rental Service

Available for NO Downtime!

**CONTACT YOUR DEALER TODAY!**

[DEALER TAG HERE]



Denterpriseintl.com