


MOBILEX[®]

Lightweight & Handheld



REVIEW ED BY:



With DOUBLElead-lined housing and lead-infused shield for maximum protection



Costs less than other brands



Take x-rays in any operator!

Save time and money with MobileX!

TECHNICAL SPECIFICATIONS

Tube Voltage	70kV
Tube Current	2mA
Input Voltage	100-240V
Frequency	50-60Hz
X-ray Tube Head	Toshiba (Japan)
Focal Spot	0.4mm
Exposure Time Range	0.04-2.0s
Weight	4.7 lbs
Battery	Li-ion
#of Exposures/charge	400



HIGH FREQUENCY
DC X-RAY
GENERATOR



USED ON SENSORS,
PHOSPHOR PLATES



LIGHTWEIGHT
4.7 LBS

OTHER FEATURES:

- Oil-cooled tube head
- Exposure time adjustment with tooth chart presets
- LED back-lit screen. Displays: time, mode, battery life, ready/indicator
- Includes neck strap and hand strap
- Can be mounted on a tripod or a pistol grip
- Rechargeable battery with 1-year warranty. Battery available to purchase and easily replaceable by user



CLEAR
IMAGES

OPTIONAL ACCESSORIES:

- Remote trigger switch w/ 6ft cord
- Pistol grip trigger switch



Outstanding Dental Equipment Since 1989

www.denterpriseintl.com | sales@denterpriseintl.com | 877.509.3180
156 East Granada Boulevard, Ormond Beach, FL 32176