

REVOLUTIONIZE VETERINARY
DENTAL IMAGING WITH

 **VET** MOBILEX



REVIEWED BY:

 **Gordon J. Christensen**
Clinicians Report[®]
A Publication of CR Foundation[®]
FEB 2019 VOL 12 ISSUE 2

 **FDA**
CLEARED

K180561

The
RELIABLE
Handheld Solution for Vets

Perfect solution for veterinary dentists seeking precision, portability, and efficiency in their practice. Its fast exposure times and adjustable settings make it ideal for a variety of dental procedures, while the durable build ensures reliable performance in demanding clinical environment.

Think Outside the Box & Join Thousands of Handheld Users!

Take x-rays in any operator!

Save time and money with MobileX!

TECHNICAL SPECIFICATIONS

Tube Voltage	70kV
Tube Current	2mA
Input Voltage	100-240V
Frequency	50-60Hz
X-ray Tube Head	Toshiba (Japan)
Focal Spot	0.4mm
Exposure Time Range	0.04-2.0s
Weight	4.7 lbs
Battery	Li-ion
#of Exposures/charge	400



HIGH FREQUENCY
DC X-RAY
GENERATOR



USED ON SENSORS,
PHOSPHOR PLATES



LIGHTWEIGHT
4.7 LBS



CLEAR
IMAGES

OTHER FEATURES:

- Oil-cooled tube head
- Exposure time adjustment with tooth chart presets
- LED back-lit screen. Displays: time, mode, battery life, ready/indicator
- Includes neck strap and hand strap
- Can be mounted on a tripod or a pistol grip
- Rechargeable battery with 1-year warranty. Battery available to purchase and easily replaceable by user

OPTIONAL ACCESSORIES:

- Remote trigger switch w/ 6ft cord
- Pistol grip trigger switch



Outstanding Dental Equipment Since 1989

www.denterpriseintl.com | sales@denterpriseintl.com | 877.509.3180
156 East Granada Boulevard, Ormond Beach, FL 32176