

# QUICKRAY<sup>DX</sup>



## The Affordable **HIGH RESOLUTION** Intraoral Sensor

**Low Cost, Never Any Hidden Fees**

The QuickRayDX Digital sensor offers superb image quality at a fraction of the cost of others. No need to sacrifice quality for price when you can have both with QuickRayDX.

COMPATIBLE WITH



QuickRay DX includes a built in TWAIN allowing for integration into all imaging software that is TWAIN compatible.





# Your Practice Deserves To Purchase The Right Sensor!

Let Our Team of Experts Help You Make the Best Choice

## SENSOR TECHNICAL DATA

SENSOR TECH	CMOS
RESOLUTION	Line Pairs >25 lp/mm
ACTIVE AREA	Size 1: 25.4 x 36.8 x 4.4 mm Size 2: 30.4 x 41.9 x 4.4 mm
SCINTILLATOR	CSI
PIXEL SIZE	1300x1800
COMPATIBILITY	Windows 32/64 Bit
DATA INTERFACE	USB 2.0 & 3.0



## QUICKRAY<sup>DX</sup> SENSOR



- Proven Durable Engineering
- Hermetically Sealed Housing
- Shock Absorbing Material
- Ergonomic Design
- Thin with Rounded Corners for Patient Comfort

## CORD



- 2M Cable Length
- Kevlar Reinforced Virtually Unbreakable
- Replaceable
- Direct High Speed 2.0 USB Connectivity
- 70.000+ bendings RA test passed



ACTUAL IMAGE

## SENSOR KIT INCLUDES



Protective Sheaths



Uni-Verse-All Positioning Kit



Wall Mount Holder



## Warranty

2 Year Warranty  
5 Year Replacement Plan



## Expedited Rental Service

Available for NO Downtime!

**CONTACT YOUR DEALER TODAY!**



[DEALER TAG HERE]