

QUICKRAY^{PRO}



NEW GENERATION Intraoral Dental X-RAY SENSOR

Unmatched Image Quality, Unbelievable Value

The **Quickray^{PRO}** digital x-ray sensor offers the best resolution on the market, backed by a market leading warranty. For high precision diagnosis, patient comfort & reliability choose **Quickray^{PRO}**!



TWAIN compatible, Quickray^{PRO} works in all leading imaging software and can be used in conjunction with existing sensors.

COMPATIBLE WITH






Your Practice Deserves To Purchase The Right Sensor! Let Our Team of Experts Help You Make the Best Choice

SENSOR TECHNICAL DATA

SENSOR TECH	CMOS
RESOLUTION	Line Pairs=35.7 lp/mm
SENSOR SIZES	Size 1: 36.2 x 24.1 x 5.9 mm Size 2: 42.8 x 30.5 x 5.7 mm
ACTIVE AREA	Size 1: 29.5 x 19.7 mm Size 2: 36.0 x 25.8 mm
PIXEL SIZE	14x14um
COMPATIBILITY	Windows 7/ 10, USB 2.0
BODY THICKNESS	5.7MM
SCINTILLATOR	CSI



QUICKRAY_{PRO} SENSOR FEATURES



ACTUAL IMAGE



CORD

- Optimized image quality at low dose
- Comfortable design
- Ease of integration
- Outstanding reliability
- Available in size 1 & 2



- 2M Cable Length
- Kevlar Reinforced Virtually Unbreakable
- Replaceable
- Direct High Speed 2.0 USB Connectivity
- 100,000+ bendings RA test passed

SENSOR KIT INCLUDES



Protective Sheaths



Uni-Verse-All Positioning Kit



Wall Mount Holder



Warranty

5 Year Warranty
7 Year Replacement Plan



Expedited Rental Service

Available for NO Downtime!

CONTACT US TODAY!



[DEALER TAG HERE]