

QUICKSCANIOS

v3 PRO

Autoclavable Tip
With spontaneous heating and windless anti-fog system

Calibration Free
Wand works with software for autocalibration

Lightweight Wand
Wand is designed for easy use. Less hand fatigue!

Secure Cable
Cable is locked in the wand for continuous operation. Never worry about cable detachment while scanning.



Your Practice Deserves To Purchase the Right Scanner!

Let Our Team of Experts Help You Make the Best Choice



QuickScan IOS:

- Scanner Size: 240 X 49.5 X 30.5mm
- Scanner Weight: 210 g
- Depth of Field: 0-15mm
- Accuracy
 - Crown: <10um
 - Full Arch: <20um
- Scan Speed
 - Single Arch: 30s
 - Full Arch with Occlusion: 90s
- Tips: Autoclavable up to 100x
- Output: STL, PLY, OBJ
- Connectivity: USB 3.0 or USB C



QuickScan IOS Software:

- AI Smart Scan
- Scan Body Matching
- Metal Scan Function
- Margin Line Drawing
- Model Builder
- Undercut Inspection
- Occlusion Distance
- Cloud Platform



Outstanding Dental Equipment Since 1989

www.denterpriseintl.com | sales@denterpriseintl.com | 877.509.3180
156 East Granada Boulevard, Ormond Beach, FL 32176