

iCAM HD⁺ plus⁺

1/3" HD CMOS
SENSOR

DUAL
CAPTURE
BUTTONS




ZOOM BUTTON

LIGHT INTENSITY
ADJUSTMENT
BUTTON




30 FRAMES PER
SECOND



WINDOWS
COMPATIBLE



1080P
16:9




CLEAR
IMAGES

Your Practice Deserves To Purchase the Right Camera!

Let Our Team of Experts Help You Make the Best Choice

TECHNICAL SPECIFICATIONS

Angle of Lens	105degrees
Camera Focus	Fixed
Depth of Field	3-60mm
Lighting Source	12 pcs LED (5600k)
Warranty	1year
Footprint	less than 6sq in
Interface	USB- UVC, TWAIN
Resolution	1080p, 720p if zoomed
Cable Length	8 feet



OTHER FEATURES:

- Workswith most major dental imaging software!
- Detachable cable
- Zoom function- up to 4x
- Adjustable LEDillumination
- Hassealed lens to prevent water absorption
- Frame rate: 30FPS @1080p
- Includes camera sleeves, wall mount holder



60 DAYS FREE
TECH SUPPORT



Outstanding Dental Equipment Since 1989

www.denterpriseintl.com | sales@denterpriseintl.com | 877.509.3180
156 East Granada Boulevard, Ormond Beach, FL 32176